

### **Protocol for making 500 mg of photo-beads**

By Huilin Zhou, updated 7/29/03

- 1) 500 mg of amino glass beads (Sigma, G4518) should have 400  $\mu\text{mol/g}$  free amino groups. Load in an empty Bio-Rad column (cat 731-1550) or whatever suitable size column. Wash once by one column volume of anhydrous DMF (Aldrich).
- 2) **Link:** Dissolve 600  $\mu\text{mol}$  each of HoBt (MW=153, in fridge, Nova Biochem, 01-62-0008) and Fmoc-aminoethyl photolinker (MW=520, in drawer, Nova Biochem, 01-60-0042) in 4 ml of dry DMF completely. Add to this solution 600  $\mu\text{mol}$  of diisopropylcarbodiimide (MW=126,  $d=0.8$ , in desiccator, Aldrich, D12540-7) for 30min. (Anytime when light-sensitive reagent is used, protect from light.)
- 3) Add 2) to the beads in 1) for 90 min incubation.
- 4) Wash beads by 3 column volumes of DMF and 2 column volumes of dry dichloromethane (Aldrich). Always remove excess solvent.
- 5) **Block:** Prepare a 2 ml mixture of 20% acetic anhydride (in desiccator, Aldrich), 30% pyridine (in hood), and 50% dichloromethane (hood). Add this mixture to the beads for 30min.
- 6) Wash beads by 2 column volumes of dichloromethane and 3 column volumes of DMF. Remove excess DMF.
- 7) **De-protect:** Make 4 ml of 20% piperidine (hood) in DMF (All common organics are from Aldrich). Incubate beads in 2 ml of 20% piperidine/DMF for 30 min. Reserve 2 ml for blank. Collect all of the 2 ml flowthrough. Use 20% piperidine/DMF as blank and measure  $A_{290}$  (absorbance at 290nm) of the 1/500 dilution (dilute by 20% piperidine/DMF) of the flowthrough. The  $A_{290}$  should be between 0.3 and 0.4. Record the value. Calculate capacity:  $(A_{290} \times \text{dilution factor} \times \text{flowthrough volume (ml)}) / (1.65 \times \text{weight of beads (mg)}) = \text{capacity } (\mu\text{mol/g})$  (should be  $\sim 400$ )
- 8) Wash beads by 5 column volumes of dry DMF again. Now the beads are ready for the next round of attachment of amino acid.
- 9) Repeat step 2 to 8 with an isotope-tagged linker other than the Fmoc-photolinker. For example, Fmoc-d0-GABA (Nova Biochem, 04-12-1088) or Fmoc-d6-GABA (protocol below).
- 10) **Attach iodoacetyl group:** Dissolve 600  $\mu\text{mol}$  of iodoacetic anhydride (MW=354, fridge, Aldrich, light sensitive) in suitable volume of dry DMF (about 4 ml) before adding to the beads. Immediately after, add in 660  $\mu\text{mol}$  of diisopropyl ethyl amine (FW=129,  $d=0.74$ , dessicator, Aldrich, D12580-6) to the beads mixture and mix it well. Let it incubate for 90 min.
- 11) Wash the beads with 5 column volumes of DMF and excess methanol. Dry in covered Bio-Spin columns in the Speedvac. Now the d0/d6-GABA-photo-beads are ready to use.

### **Protocol for Fmoc-d6-GABA**

(600  $\mu\text{mol}$  for 500 mg beads) (Andy Tao's improved protocol)

- 1) Dissolve 600  $\mu\text{mol}$  d6-GABA (MW=109, shelf, Isotec) in 3 ml 10% sodium carbonate in  $\text{H}_2\text{O}$ , in glass vial with stir bar.
- 2) Dissolve 600  $\mu\text{mol}$  Fmoc-Osu (MW=337, -20C, Nova Biochem) in 3 ml dimethoxyethane. Add drop-wise to d6-GABA in vial. Stir at least 2 hours.
- 3) Add 20 ml  $\text{H}_2\text{O}$ , making clear solution. Centrifuge to remove precipitate.
- 4) Add up to 600  $\mu\text{l}$  conc. HCl to pH  $\sim 2$ . White precipitate forms.

- 5) Extract Fmoc-GABA with 3-4 ml ethyl acetate; collect organic phase and retain aqueous to repeat extraction 4-5 times, combining extracts.
- 6) Wash extract with H<sub>2</sub>O.
- 7) Weigh 4 empty microtubes and record weights. Dry extract (about 15 ml) in microtubes, adding more to tubes as it dries. Record final weights. Compare to calculated weight of 600  $\mu$ mol of d6-Fmoc-GABA (MW=331) = 199 mg.

Step	beads reagent	100 mg 120 $\mu$ mol	500 mg 600 $\mu$ mol
<b>Attach photolinker</b>	DMF	800 $\mu$ l	4 ml
	HoBt	18.4 mg	92 mg
	Fmoc photo-linker	62.4 mg	312 mg
	DIC	18.9 $\mu$ l	95 $\mu$ l
<b>Block</b>	Dichloromethane	500 $\mu$ l	1 ml
	Acetic anhydride	200 $\mu$ l	400 $\mu$ l
	Pyridine	300 $\mu$ l	600 $\mu$ l
<b>De-protect</b>	DMF	1.6 ml	3.2 ml
	Piperidine	400 $\mu$ l	800 $\mu$ l
<b>Link GABA</b>	DMF	800 $\mu$ l	4 ml
	HoBt	18.4 mg	92 mg
	Fmoc-d0-GABA	39 mg	195 mg
	DIC	18.9 $\mu$ l	95 $\mu$ l
<b>Attach iodo</b>	DMF	800 $\mu$ l	4 ml
	Iodoacetic anhydride	42.5 mg	212 mg
	Diisopropyl ethyl amine	23 $\mu$ l (132 $\mu$ mol)	115 $\mu$ l (660 $\mu$ mol)
<b>Fmoc-GABA</b>	Sodium carbonate/H <sub>2</sub> O		300 mg/3 ml
	d6-GABA		65 mg
	Fmoc-Osu/DMF		202 mg/3 ml